


Cabinet Assembly/Dimensions




Clear Inside Drawer Dimensions:		Clear Inside Cupboard Dimensions	
Typical Depth	18-1/8" (Front to Back)	Typical Depth	20-7/8" (Front to Back)
3-5/8"H Drawer	2-7/8" High	Full Height	29"
7-1/4"H Drawer	6-1/2" High	2 Drawer Height	14-1/2"
14-1/2"H Drawer	13-3/4" High	3 Drawer Height	21-3/4"
18"W Cabinet	14-3/4" Wide	18"W Cabinet	14-5/8" Total Inside; 17-7/8" Clear Inside
24"W Cabinet	20-3/4" Wide	24"W Cabinet	20-5/8" Total Inside; 23-7/8" Clear Inside
35"W Cabinet	31-3/4" Wide for full width drawer or (2) at 14-3/4"	35"W Cabinet	31-5/8" Total Inside (Upright Removed); 34-7/8" Clear Inside
47"W Cabinet	43-3/4" Wide for full width drawer or (2) at 20-3/4"	47"W Cabinet	43-5/8" Total Inside (Upright Removed); 46-7/8" Clear Inside
			Single Door Combo Unit 14-5/8" Total Inside
			Single Door Combo Unit 20-5/8" Total Inside

Standing Height Base Cabinets

7000 Series - Single Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:

- * Specify swing for cupboard doors on *single* units by substituting '-' with 'R' or 'L' after series number (Units are shown with *left* swing). Ex. 7009R24 for right swing cupboard unit, 24" wide.
- * For finished backs, add suffix 'FB' to end of model number (Ex. 7001R24FB).
- * Specify locks if required.




Model 7001

Unit includes full height cupboard with one adjustable shelf.

7001-18 = 18" Wide

7001-24 = 24" Wide



Model 7005

Unit includes two 3-5/8" drawers above three 7-1/4" drawers.

7005-18 = 18" Wide

7005-24 = 24" Wide




Model 7002

Unit includes one 7-1/4" drawer above cupboard with one adjustable shelf.

7002-18 = 18" Wide

7002-24 = 24" Wide




Model 7006

Unit includes eight 3-5/8" drawers.

7006-18 = 18" Wide

7006-24 = 24" Wide




Model 7004

Unit includes four 7-1/4" drawers.

7004-18 = 18" Wide

7004-24 = 24" Wide



Model 7008

Open unit includes one adjustable shelf.

7008-18 = 18" Wide


7008-24 = 24" Wide

Standing Height Base Cabinets

7000 Series - Single Units

Dimensions:

End View



Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.
- * All file drawers have follow-blocks and full extension slides.

Ordering Information:

- * Specify swing for cupboard doors on *single* units by substituting '-' with 'R' or 'L' after series number (Units are shown with *left* swing). Ex. 7009R24 for right swing cupboard unit, 24" wide.
- * For finished backs, add suffix 'FB' to end of model number (Ex. 7010-24FB).
- * Specify locks if required.




Model 7009

Unit includes two 3-5/8" drawers above cupboard with one adjustable shelf.

7009-18 = 18" Wide

7009-24 = 24" Wide



Model 7013

Unit includes two 14-1/2" file drawers.

7013-18 = 18" Wide

7013-24 = 24" Wide




Model 7010

Unit includes four 3-5/8" drawers above two 7-1/4" drawers.

7010-18 = 18" Wide

7010-24 = 24" Wide




Model 7014

Unit includes four 3-5/8" drawers above 14-1/2" file drawer.

7014-18 = 18" Wide

7014-24 = 24" Wide




Model 7011

Unit includes two 7-1/4" drawers above 14-1/2" file drawer.

7011-18 = 18" Wide

7011-24 = 24" Wide




Model 7012

Unit includes two 3-5/8" drawers above 7-1/4" drawer and 14-1/2" file drawer.

7012-18 = 18" Wide

7012-24 = 24" Wide



Model 7017

Trash unit includes tilt-out door with attached waste container.

7017-18 = 18" Wide


7017-24 = 24" Wide

Standing Height Base Cabinets

7100 Series - Double & Combination Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:

- * For finished backs, add suffix 'FB' to model number.
- * Specify locks if required.



Model 7101

Unit includes one full width adjustable shelf.


- | | |
|--------------------|--------------------|
| 7101-29 = 29" Wide | 7101-41 = 41" Wide |
| 7101-35 = 35" Wide | 7101-47 = 47" Wide |



Model 7102

Unit includes two 7-1/4" drawers above two cupboards, and one full width adjustable shelf.

- | | |
|--------------------|--------------------|
| 7102-29 = 29" Wide | 7102-41 = 41" Wide |
| 7102-35 = 35" Wide | 7102-47 = 47" Wide |



Model 7104

Unit includes eight 7-1/4" drawers.

- | | |
|--------------------|--------------------|
| 7104-29 = 29" Wide | 7104-41 = 41" Wide |
| 7104-35 = 35" Wide | 7104-47 = 47" Wide |



Model 7105

Unit includes four 3-5/8" drawers above six 7-1/4" drawers.


- | | |
|--------------------|--------------------|
| 7105-29 = 29" Wide | 7105-41 = 41" Wide |
| 7105-35 = 35" Wide | 7105-47 = 47" Wide |

Standing Height Base Cabinets

7100 Series - Double & Combination Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:


- * Specify right or left side cupboard door on *combination* units by substituting '-' with 'L' or 'R' after series number (Units are shown with right side doors) Ex. 7110L35 for left side door.
- * For finished backs, add suffix 'FB' to model number.
- * Specify locks if required.



Model 7106

Unit includes sixteen 3-5/8" drawers.


7106-29 = 29" Wide	7106-41 = 41" Wide
7106-35 = 35" Wide	7106-47 = 47" Wide



Model 7108

Unit includes one full width adjustable shelf.


7108-29 = 29" Wide	7108-41 = 41" Wide
7108-35 = 35" Wide	7108-47 = 47" Wide



Model 7109

Unit includes four 3-5/8" drawers above two cupboards, and one full width adjustable shelf.

7109-29 = 29" Wide	7109-41 = 41" Wide
7109-35 = 35" Wide	7109-47 = 47" Wide



Model 7110

Combination unit includes four 3-5/8" drawers above two 7-1/4" drawers, and full height cupboard with one half-width adjustable shelf.


7110-29 = 29" Wide	7110-41 = 41" Wide
7110-35 = 35" Wide	7110-47 = 47" Wide

Standing Height Base Cabinets

7100 Series - Double & Combination Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.
- * All file drawers have follow-blocks and full extension slides.

Ordering Information:

- * Specify right or left side cupboard door on *combination* units by substituting '-' with 'L' or 'R' after series number (Units are shown with right side doors) Ex. 7111L35 for left side door.
- * For finished backs, add suffix 'FB' to model number.
- * If locks are required, specify quantity and location (drawer or door) for each unit.



Model 7111


Combination unit includes two 7-1/4" drawers above 14-1/2" file drawer, and full height cupboard with one half-width adjustable shelf.

7111-29 = 29" Wide

7111-41 = 41" Wide

7111-35 = 35" Wide

7111-47 = 47" Wide



Model 7113


Combination unit includes four 14-1/2" file drawers.

7113-29 = 29" Wide

7113-41 = 41" Wide

7113-35 = 35" Wide

7113-47 = 47" Wide



Model 7114


Combination unit includes four 3-5/8" drawers above one 14-1/2" file drawer, and full height cupboard with one half-width adjustable shelf.

7114-29 = 29" Wide

7114-41 = 41" Wide

7114-35 = 35" Wide

7114-47 = 47" Wide



Model 7119

Unit includes one 7-1/4" full-width drawer above two cupboards, and one full width adjustable shelf.

7119-29 = 29" Wide

7119-41 = 41" Wide

7119-35 = 35" Wide


7119-47 = 47" Wide

Standing Height Base Cabinets

7100 Series - Double & Combination Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:

- * Specify right or left side cupboard door on *combination* units by substituting '-' with 'L' or 'R' after series number (Units are shown with right side doors) Ex. 7121L35 for left side door.
- * For finished backs, add suffix 'FB' to model number.
- * If locks are required, specify quantity and location (drawer or door) for each unit.



Model 7121


Combination unit includes four 7-1/4" drawers, and full height cupboard with one half-width adjustable shelf.

7121-29 = 29" Wide

7121-41 = 41" Wide

7121-35 = 35" Wide

7121-47 = 47" Wide



Model 7122


Combination unit includes two 3-5/8" drawers above three 7-1/4" drawers, and full height cupboard with one half-width adjustable shelf.

7122-29 = 29" Wide

7122-41 = 41" Wide

7122-35 = 35" Wide

7122-47 = 47" Wide



Model 7124

Combination unit includes eight 3-5/8" drawers, and full height cupboard with one half-width adjustable shelf.

7124-29 = 29" Wide

7124-41 = 41" Wide

7124-35 = 35" Wide

7124-47 = 47" Wide



Model 7125

Combination unit includes four 7-1/4" drawers, and two 3-5/8" drawers above cupboard with one half-width adjustable shelf.

7125-29 = 29" Wide

7125-41 = 41" Wide

7125-35 = 35" Wide


7125-47 = 47" Wide

Standing Height Base Cabinets

7100 Series - Double & Combination Units

Dimensions:

End View



Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:

- * Specify right or left side cupboard door on *combination* units by substituting '-' with 'L' or 'R' after series number (Units are shown with right side doors) Ex. 7126L35 for left side door.
- * For finished backs, add suffix 'FB' to model number.
- * If locks are required, specify quantity and location (drawer or door) for each unit.



Model 7126

Combination unit includes five 7-1/4" drawers, with one drawer above cupboard with one half-width adjustable shelf.

7126-29 = 29" Wide

7126-41 = 41" Wide

7126-35 = 35" Wide

7126-47 = 47" Wide



Model 7127


Combination unit includes ten 3-5/8" drawers, with two drawers above cupboard with one half-width adjustable shelf.

7127-29 = 29" Wide

7127-41 = 41" Wide

7127-35 = 35" Wide

7127-47 = 47" Wide



Model 7130

Unit includes four 7-1/4" full-width drawers.

7130-29 = 29" Wide

7130-41 = 41" Wide

7130-35 = 35" Wide

7130-47 = 47" Wide



Model 7131

Unit includes two 7-1/4" drawers above three 7-1/4" full width drawers.

7131-29 = 29" Wide

7131-41 = 41" Wide

7131-35 = 35" Wide


7131-47 = 47" Wide

Standing Height Base Cabinets

7100 Series - Double & Combination Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:

- * For finished backs, add suffix 'FB' to model number.
- * If locks are required, specify quantity and location (drawer or door) for each unit.



Model 7132


Unit includes four 7-1/4" drawers above two 7-1/4" full width drawers.

7132-29 = 29" Wide

7132-41 = 41" Wide

7132-35 = 35" Wide

7132-47 = 47" Wide



Model 7136


Dead corner unit includes full height cupboard and fixed panel, and one full width adjustable shelf. Specify **Left** or **Right** side cupboard. (Shown with Left Side Cupboard)

7136-29 = 29" Wide

7136-41 = 41" Wide

7136-35 = 35" Wide

7136-47 = 47" Wide



Model 7141


Unit includes full height sliding door cupboard and one full width adjustable shelf.

7141-29 = 29" Wide

7141-41 = 41" Wide

7141-35 = 35" Wide

7141-47 = 47" Wide



Model 7143

Unit includes one 7-1/4" full-width drawer above sliding door cupboard, and one full width adjustable shelf.

7143-29 = 29" Wide

7143-41 = 41" Wide

7143-35 = 35" Wide


7143-47 = 47" Wide

Standing Height Base Cabinets

7100 Series - Double & Combination Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.
- * All file drawers have follow-blocks and full extension slides.

Ordering Information:

- * Specify right or left side cupboard door on *combination* units by substituting '-' with 'L' or 'R' after series number (Units are shown with right side doors) Ex. 7149L35 for left side door.
- * For finished backs, add suffix 'FB' to model number.
- * If locks are required, specify quantity and location (drawer or door) for each unit.



Model 7145

Unit includes two 7-1/4" drawers above sliding door cupboard, and one full width adjustable shelf.

7145-29 = 29" Wide

7145-41 = 41" Wide

7145-35 = 35" Wide

7145-47 = 47" Wide



Model 7146


Combination unit includes four 3-5/8" drawers above sliding door cupboard, and full width adjustable shelf.

7146-29 = 29" Wide

7146-41 = 41" Wide

7146-35 = 35" Wide

7146-47 = 47" Wide



Model 7149


Combination unit includes two 3-5/8" drawers above 7-1/4" h drawer and 14-1/2" file drawer with follow-block, and full height cupboard with one half-width adjustable shelf.

7149-29 = 29" Wide

7149-41 = 41" Wide

7149-35 = 35" Wide

7149-47 = 47" Wide



Model 7167

Combination unit includes two 14-1/2" file drawers with follow-blocks, and full height cupboard with one half-width adjustable shelf.

7167-29 = 29" Wide

7167-41 = 41" Wide

7167-35 = 35" Wide


7167-47 = 47" Wide

Sitting Height Base Cabinets

7300 Series - Single Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:

Specify swing for cupboard doors on *single* units by substituting '-' with 'R' or 'L' after series number (Units are shown with *left* swing). Ex. 7301R24 for right swing cupboard unit, 24" wide).


- * For finished backs, add suffix 'FB' to model number.
- * Specify locks if required.



Model 7301

Unit includes full height cupboard with one adjustable shelf.


- 7301-18 = 18" Wide
- 7301-24 = 24" Wide



Model 7305

Unit includes two 3-5/8" drawers above two 7-1/4" drawers.

- 7305-18 = 18" Wide
- 7305-24 = 24" Wide



Model 7302

Unit includes one 7-1/4" drawer above cupboard with one adjustable shelf.


- 7302-18 = 18" Wide
- 7302-24 = 24" Wide



Model 7306

Unit includes six 3-5/8" drawers.


- 7306-18 = 18" Wide
- 7306-24 = 24" Wide



Model 7304

Unit includes three 7-1/4" drawers.

- 7304-18 = 18" Wide
- 7304-24 = 24" Wide



Model 7308

Open unit includes one adjustable shelf.


- 7308-18 = 18" Wide
- 7308-24 = 24" Wide

Sitting Height Base Cabinets

7300 Series - Single Units

Dimensions:

End View



Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.
- * All file drawers have follow-blocks and full extension slides.

Ordering Information:

Specify swing for cupboard doors on *single* units by substituting '-' with 'R' or 'L' after series number (Units are shown with *left* swing). Ex. 7309**R**24 for right swing cupboard unit, 24" wide


- * For finished backs, add suffix 'FB' to model number.
- * Specify locks if required.



Model 7309

Unit includes two 3-5/8" drawers above cupboard with one adjustable shelf.

- 7309-18 = 18" Wide
- 7309-24 = 24" Wide



Model 7313

Unit includes two 10-7/8" file drawers.


- 7313-18 = 18" Wide
- 7313-24 = 24" Wide



Model 7310

Unit includes four 3-5/8" drawers above one 7-1/4" drawer.

- 7310-18 = 18" Wide
- 7310-24 = 24" Wide



Model 7314

Unit includes two 3-5/8" drawers above 14-1/2" file drawer.

- 7314-18 = 18" Wide
- 7314-24 = 24" Wide



Model 7311

Unit includes one 7-1/4" drawer above 14-1/2" file drawer.


- 7311-18 = 18" Wide
- 7311-24 = 24" Wide

Sitting Height Base Cabinets

7400 Series - Double & Combination Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:

- * For finished backs, add suffix 'FB' to model number.
- * Specify locks if required.



Model 7401


Unit includes one full width adjustable shelf.

7401-29 = 29" Wide

7401-41 = 41" Wide

7401-35 = 35" Wide

7401-47 = 47" Wide



Model 7402

Unit includes two 7-1/4" drawers above two cupboards.

7402-29 = 29" Wide

7402-41 = 41" Wide

7402-35 = 35" Wide

7402-47 = 47" Wide



Model 7404

Unit includes six 7-1/4" drawers.

7404-29 = 29" Wide

7404-41 = 41" Wide

7404-35 = 35" Wide

7404-47 = 47" Wide



Model 7405


Unit includes four 3-5/8" drawers above four 7-1/4" drawers.

7405-29 = 29" Wide

7405-41 = 41" Wide

7405-35 = 35" Wide

7405-47 = 47" Wide



Model 7406

Unit includes twelve 3-5/8" drawers.

7406-29 = 29" Wide

7406-41 = 41" Wide

7406-35 = 35" Wide


7406-47 = 47" Wide

Sitting Height Base Cabinets

7400 Series - Double & Combination Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.
- * All file drawers have follow-blocks and full extension slides.

Ordering Information:

- * Specify right or left side cupboard door on *combination* units by substituting '-' with 'L' or 'R' after series number (Units are shown with right side doors) Ex. 7411L35 for left side door.
- * For finished backs, add suffix 'FB' to model number.
- * Specify locks if required.



Model 7408

Unit includes one full width adjustable shelves.

- | | |
|--------------------|--------------------|
| 7408-29 = 29" Wide | 7408-41 = 41" Wide |
| 7408-35 = 35" Wide | 7408-47 = 47" Wide |



Model 7409

Unit includes four 3-5/8" drawers above two cupboards.


- | | |
|--------------------|--------------------|
| 7409-29 = 29" Wide | 7409-41 = 41" Wide |
| 7409-35 = 35" Wide | 7409-47 = 47" Wide |



Model 7410

Combination unit includes four 3-5/8" drawers above one 7-1/4" drawers, and full height cupboard with one half-width adjustable shelf.


- | | |
|--------------------|--------------------|
| 7410-29 = 29" Wide | 7410-41 = 41" Wide |
| 7410-35 = 35" Wide | 7410-47 = 47" Wide |



Model 7411

Combination unit includes one 7-1/4" drawer above 14-1/2" file drawer, and full height cupboard with one half-width adjustable shelf.

- | | |
|--------------------|--------------------|
| 7411-29 = 29" Wide | 7411-41 = 41" Wide |
| 7411-35 = 35" Wide | 7411-47 = 47" Wide |



Model 7414

Combination unit includes two 3-5/8" drawers above one 14-1/2" file drawer, and full height cupboard with one half-width adjustable shelf.


- | | |
|--------------------|--------------------|
| 7414-29 = 29" Wide | 7414-41 = 41" Wide |
| 7414-35 = 35" Wide | 7414-47 = 47" Wide |

Sitting Height Base Cabinets

7400 Series - Double & Combination Units

Dimensions:

End View



Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:

- * Specify right or left side cupboard door on *combination* units by substituting '-' with 'L' or 'R' after series number (Units are shown with right side doors) Ex. 7425L35 for left side door.
- * For finished backs, add suffix 'FB' to model number.
- * Specify locks if required



Model 7419

Unit includes one 7-1/4" full-width drawer above two cupboards.

7419-29 = 29" Wide

7419-41 = 41" Wide

7419-35 = 35" Wide

7419-47 = 47" Wide



Model 7421

Combination unit includes three 7-1/4" drawers, and full height cupboard with one half-width adjustable shelf. (Unit shown with left side door)

7421-29 = 29" Wide

7421-41 = 41" Wide

7421-35 = 35" Wide

7421-47 = 47" Wide



Model 7422

Combination unit includes two 3-5/8" drawers above two 7-1/4" drawers, and full height cupboard with one half-width adjustable shelf.

7422-29 = 29" Wide

7422-41 = 41" Wide

7422-35 = 35" Wide

7422-47 = 47" Wide



Model 7424

Combination unit includes six 3-5/8" drawers, and full height cupboard with one half-width adjustable shelf.

7424-29 = 29" Wide

7424-41 = 41" Wide

7424-35 = 35" Wide

7424-47 = 47" Wide



Model 7425

Combination unit includes three 7-1/4" drawers, and two 3-5/8" drawers above cupboard with one half-width adjustable shelf.

7425-29 = 29" Wide

7425-41 = 41" Wide

7425-35 = 35" Wide


7425-47 = 47" Wide

Sitting Height Base Cabinets

7400 Series - Double & Combination Units

Dimensions:

End View



Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:

- * Specify right or left side cupboard door on *combination* units by substituting '-' with 'L' or 'R' after series number (Units are shown with right side doors) Ex. 7426**L**35 for left side door.
- * For finished backs, add suffix '**FB**' to model number.
- * Specify locks if required.



Model 7426

Combination unit includes four 7-1/4" drawers, with one drawer above cupboard with one half-width adjustable shelf.

7426-29 = 29" Wide	7426-41 = 41" Wide
7426-35 = 35" Wide	7426-47 = 47" Wide



Model 7427

Combination unit includes eight 3-5/8" drawers, with two drawers above cupboard with one half-width adjustable shelf.

7427-29 = 29" Wide	7427-41 = 41" Wide
7427-35 = 35" Wide	7427-47 = 47" Wide



Model 7430

Unit includes three 7-1/4" full-width drawers.


7430-29 = 29" Wide	7430-41 = 41" Wide
7430-35 = 35" Wide	7430-47 = 47" Wide



Model 7431

Unit includes two 7-1/4" drawers above two 7-1/4" full width drawers.

7431-29 = 29" Wide	7431-41 = 41" Wide
7431-35 = 35" Wide	7431-47 = 47" Wide



Model 7432

Unit includes four 7-1/4" drawers above one 7-1/4" full width drawer.


7432-29 = 29" Wide	7432-41 = 41" Wide
7432-35 = 35" Wide	7432-47 = 47" Wide

Sitting Height Base Cabinets

7400 Series - Double & Combination Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:

- * For finished backs, add suffix '**FB**' to model number.
- * Specify locks if required.



Model 7436

Dead corner unit includes full height cupboard and fixed panel, and one full width adjustable shelf. Specify **Left** or **Right** cupboard. (Shown with Left Side cupboard)

7436-29 = 29" Wide

7436-41 = 41" Wide

7436-35 = 35" Wide

7436-47 = 47" Wide



Model 7443

Unit includes one 7-1/4" full-width drawer above sliding door cupboard.

7443-29 = 29" Wide

7443-41 = 41" Wide

7443-35 = 35" Wide

7443-47 = 47" Wide



Model 7445

Unit includes two 7-1/4" drawers above sliding door cupboard.

7445-29 = 29" Wide

7445-41 = 41" Wide

7445-35 = 35" Wide

7445-47 = 47" Wide



Model 7446

Combination unit includes four 3-5/8" drawers above sliding door cupboard.

7446-29 = 29" Wide

7446-41 = 41" Wide

7446-35 = 35" Wide


7446-47 = 47" Wide

ADA Height Base Cabinets

7000 Series - Single Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:

- * Specify swing for cupboard doors on *single* units by substituting '-' with 'R' or 'L' after series number (Units are shown with left swing). Ex. 7001ADAR24 for right swing door unit, 24" wide.
- * For finished backs, add suffix 'FB' to end of model number (Ex. 7001ADAR24FB).
- * Specify locks if required.



Model 7001ADA

Unit includes full height cupboard with one adjustable shelf.

7001ADA-18 = 18" W

7001ADA-24 = 24" W




Model 7005ADA

Unit includes three 3-5/8" drawers above two 7-1/4" drawers.

7005ADA-18 = 18" W

7005ADA-24 = 24" W




Model 7002ADA

Unit includes one 7-1/4" drawer above cupboard with one adjustable shelf.

7002ADA-18 = 18" W

7002ADA-24 = 24" W




Model 7006ADA

Unit includes seven 3-5/8" drawers.

7006ADA-18 = 18" W

7006ADA-24 = 24" W




Model 7004ADA

Unit includes one 3-5/8" drawers above three 7-1/4" drawers.

7004ADA-18 = 18" W

7004ADA-24 = 24" W



Model 7008ADA

Open unit includes one adjustable shelf.

7008ADA-18 = 18" W


7008ADA-24 = 24" W

ADA Height Base Cabinets

7000 Series - Single Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.
- * All file drawers have follow-blocks and full extensions slides.

Ordering Information:


- * Specify swing for cupboard doors on *single* units by substituting '-' with 'R' or 'L' after series number (Units are shown with *left* swing). Ex. 7009ADAR24 for right swing door unit, 24" wide.
- * For finished backs, add suffix 'FB' to end of model number.
- * Specify locks if required.



Model 7009ADA

Unit includes two 3-5/8" drawers above cupboard with one adjustable shelf.


7009ADA-18 = 18" Wide
7009ADA-24 = 24" Wide



Model 7012ADA

Unit includes three 3-5/8" drawers above a 14-1/2" file drawer.


7012ADA-18 = 18" Wide
7012ADA-24 = 24" Wide



Model 7010ADA

Unit includes five 3-5/8" drawers above one 7-1/4" drawers.


7010ADA-18 = 18" Wide
7010ADA-24 = 24" Wide



Model 7013ADA

Unit includes one 10-7/8" file drawer above a 14-1/2" file drawer.


7013ADA-18 = 18" Wide
7013ADA-24 = 24" Wide



Model 7011ADA

Unit includes one 3-5/8" drawer above one 7-1/4" drawer and 14-1/2" file drawer.

7011ADA-18 = 18" Wide
7011ADA-24 = 24" Wide



Model 7014ADA

Open unit includes four 3-5/8" drawers above one 10-7/8" file drawer.


7014ADA-18 = 18" Wide
7014ADA-24 = 24" Wide

ADA Height Base Cabinets

7100 Series - Double & Combination Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.

Ordering Information:


- * Specify right or left side cupboard door by substituting '-L' with 'L' or '-R' after series number (Units are shown with right side doors) Ex. 7110ADAL35 for left side cupboard.
- * For finished backs, add suffix 'FB' to end of model number (Ex. 7110ADAR24FB).
- * Specify locks if required.



Model 7101ADA

Unit includes one full width adjustable shelf.


- 7101ADA-29 = 29" W
- 7101ADA-41 = 41" W
- 7101ADA-35 = 35" W
- 7101ADA-47 = 47" W



Model 7106ADA

Unit includes fourteen 3-5/8" drawers.


- 7106ADA-29 = 29" W
- 7106ADA-41 = 41" W
- 7106ADA-35 = 35" W
- 7106ADA-47 = 47" W



Model 7102ADA

Unit includes two 7-1/4" drawers above two cupboards, and one full width adjustable shelf.

- 7102ADA-29 = 29" W
- 7102ADA-41 = 41" W
- 7102ADA-35 = 35" W
- 7102ADA-47 = 47" W



Model 7108ADA

Unit includes one full width adjustable shelf.

- 7108ADA-29 = 29" W
- 7108ADA-41 = 41" W
- 7108ADA-35 = 35" W
- 7108ADA-47 = 47" W



Model 7104ADA

Unit includes two 3-5/8" drawers above six 7-1/4" drawers.

- 7104ADA-29 = 29" W
- 7104ADA-41 = 41" W
- 7104ADA-35 = 35" W
- 7104ADA-47 = 47" W



Model 7109ADA

Unit includes four 3-5/8" drawers above two cupboards, and one full width adjustable shelf.


- 7109ADA-29 = 29" W
- 7109ADA-41 = 41" W
- 7109ADA-35 = 35" W
- 7109ADA-47 = 47" W



Model 7105ADA

Unit includes six 3-5/8" drawers above four 7-1/4" drawers.

- 7105ADA-29 = 29" W
- 7105ADA-41 = 41" W
- 7105ADA-35 = 35" W
- 7105ADA-47 = 47" W



Model 7110ADA

Unit includes five 3-5/8" drawers above one 7-1/4" drawer, full height cupboard, one half width adjustable shelf.


- 7110ADA-29 = 29" W
- 7110ADA-41 = 41" W
- 7110ADA-35 = 35" W
- 7110ADA-47 = 47" W

ADA Height Base Cabinets

7100 Series - Double & Combination Units

Dimensions:

End View




Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.
- * All file drawers have follow-blocks and full extension slides.

Ordering Information:

- * Specify right or left side cupboard door by substituting '-' with 'L' or 'R' after series number (Units are shown with right side doors) Ex. 7111ADAL35 for left side cupboard.
- * For finished backs, add suffix 'FB' to end of model number.
- * Specify locks if required.


Model 7111ADA



Unit includes one 3-5/8" drawer above one 7-1/4" drawer, and 14-1/2" file drawer, full height cupboard, one half width adjustable shelf.

- 7111ADA-29 = 29" W
- 7111ADA-41 = 41" W
- 7111ADA-35 = 35" W
- 7111ADA-47 = 47" W


Model 7119ADA



Unit includes one 7-1/4" full-width drawer above two cupboards, one full width adjustable shelf.

- 7119ADA-29 = 29" W
- 7119ADA-41 = 41" W
- 7119ADA-35 = 35" W
- 7119ADA-47 = 47" W


Model 7113ADA



Unit includes two 10-7/8" file drawers, above two 14-1/2" file drawers.

- 7113ADA-29 = 29" W
- 7113ADA-41 = 41" W
- 7113ADA-35 = 35" W
- 7113ADA-47 = 47" W


Model 7121ADA



Unit includes one 3-5/8" drawer above three 7-1/4" drawers, full height cupboard, one half width adjustable shelf.

- 7121ADA-29 = 29" W
- 7121ADA-41 = 41" W
- 7121ADA-35 = 35" W
- 7121ADA-47 = 47" W

Model 7114ADA



Unit includes four 3-5/8" drawers above one 14-1/2" file drawer, full height cupboard, one half width adjustable shelf.

- 7114ADA-29 = 29" W
- 7114ADA-41 = 41" W
- 7114ADA-35 = 35" W
- 7114ADA-47 = 47" W

Model 7122ADA



Unit includes three 3-5/8" drawers above two 7-1/4" drawers, full height cupboard, one half width adjustable shelf.


- 7122ADA-29 = 29" W
- 7122ADA-41 = 41" W
- 7122ADA-35 = 35" W
- 7122ADA-47 = 47" W

ADA Height Base Cabinets

7100 Series - Double & Combination Units

Dimensions:

End View



Product Features:

- * Shelves are adjustable in 1/2" increments.
- * All cabinets have removable rear access panels.
- * All file drawers have follow-blocks and full extension slides.

Ordering Information:

- * Specify right or left side cupboard door by substituting '-' with 'L' or 'R' after series number (Units are shown with right side doors) Ex. 7125ADAL35 for left side cupboard.
- * For finished backs, add suffix 'FB' to end of model number.
- * Specify locks if required.


Model 7125ADA



Unit includes three 3-5/8" drawers above three 7-1/4" drawers, cupboard, one half width adjustable shelf.

- 7125ADA-29 = 29" W
- 7125ADA-41 = 41" W
- 7125ADA-35 = 35" W
- 7125ADA-47 = 47" W

Model 7143ADA



Unit includes one 7-1/4" full width drawer above sliding door cupboard, one full width adjustable shelf.

- 7143ADA-29 = 29" W
- 7143ADA-41 = 41" W
- 7143ADA-35 = 35" W
- 7143ADA-47 = 47" W


Model 7126ADA



Unit includes one 3-5/8" drawer, four 7-1/4" drawers, cupboard, one half width adjustable shelf.

- 7126ADA-29 = 29" W
- 7126ADA-41 = 41" W
- 7126ADA-35 = 35" W
- 7126ADA-47 = 47" W


Model 7145ADA



Unit includes two 7-1/4" drawers above sliding door cupboard, one full width adjustable shelf.

- 7145ADA-29 = 29" W
- 7145ADA-41 = 41" W
- 7145ADA-35 = 35" W
- 7145ADA-47 = 47" W


Model 7130ADA



Unit includes one 3-5/8" full width drawer above three 7-1/4" drawers.

- 7130ADA-29 = 29" W
- 7130ADA-41 = 41" W
- 7130ADA-35 = 35" W
- 7130ADA-47 = 47" W

Model 7146ADA



Unit includes four 3-5/8" drawers, above sliding door cupboard, one full width adjustable shelf.

- 7146ADA-29 = 29" W
- 7146ADA-41 = 41" W
- 7146ADA-35 = 35" W
- 7146ADA-47 = 47" W


Model 7131ADA



Unit includes two 7-1/4" full-drawers above one 3-5/8" full width drawer and two 7-1/4" full width drawers.

- 7131ADA-29 = 29" W
- 7131ADA-41 = 41" W
- 7131ADA-35 = 35" W
- 7131ADA-47 = 47" W

Model 7167ADA



Unit includes one 10-7/8" file drawer above a 14-1/2" file drawer, full height cupboard, one half width adjustable shelf.


- 7167ADA-29 = 29" W
- 7167ADA-41 = 41" W
- 7167ADA-35 = 35" W
- 7167ADA-47 = 47" W

Sink Base Cabinets

7000, 7100, 7300, 7400 Series

Dimensions:

End Views




Product Features:


- * Louvered front panels.
- * Low back for easy plumbing access (18" A.F.F.)
- * ADA Height Units are also available, 32-1/8"H

Ordering Information:


- * Specify swing for cupboard doors on *single* units by substituting '-' with 'R' or 'L' after series number (Unit shown with *left* swing door. Ex. 7003R18 for right swing cupboard unit, 18"W).
- * Specify sink hangers if required.
- * Specify locks if required.




Model 7003
 Standing height single sink base unit includes a louvered panel above the cupboard.
 7003-18 = 18" Wide
 7003-24 = 24" Wide




Model 7303
 Sitting height single sink base unit includes a louvered panel above the cupboard.
 7303-18 = 18" Wide
 7303-24 = 24" Wide



Model 7103
 Standing height double sink base unit includes louvered panels above the cupboards.
 7103-29 = 29" Wide
 7103-35 = 35" Wide
 7103-41 = 41" Wide
 7103-47 = 47" Wide



Model 7403
 Sitting height double sink base unit includes louvered panels above the cupboards.
 7403-29 = 29" Wide
 7403-35 = 35" Wide
 7403-41 = 41" Wide
 7403-47 = 47" Wide



End Unit Models 7103, 7403
 End sink base unit includes louvered panels above the spaced cupboard doors, and a 3-sided toe space base for end-of-island or peninsula applications.

<u>Standing Height</u>	<u>Sitting Height</u>
7103-52 = 52" Wide	7403-52 = 52" Wide
7103-58 = 58" Wide	7403-58 = 58" Wide

Corner Base Cabinets

7100, 7400 Series

Product Features:




- * Louvered front panel units 7173 and 7473.
- * Removable rear access panels.
- * Space for optional Lazy Susan shelves (24" diameter maximum) in 7171 and 7471 (see pg. 35).

Ordering Information:

- * Cupboard doors will be furnished *left-hinged* unless otherwise specified by substituting 'L' with 'R' after series number (Units shown with left swing doors). Ex. 7171**R**18 for right swing.
- * Specify sink hangers if required (for undermount sinks).
- * Specify locks if required.

Dimensions:

End Views



Models 7171/7471


Corner unit includes a full height cupboard and one adjustable shelf.

Standing Height 35-3/4"

7171-18

Sitting Height 28-1/2"

7471-18



Models 7173/7473

Corner sink base unit includes a louvered panel above the cupboard. (Shown with optional sink hangers)

Standing Height 35-3/4"

7173-18

Sitting Height 28-1/2"


7473-18

Knee Space Drawers

7200 Series – Table/Knee Space Components


Dimensions:

End View



Ordering Information:

- * Specify locks if required.
- * Specify width of each drawer End in inches (Example: 7200-18 with 4" Left End and 8" Right End).
- * If adding Legs, the minimum drawer END size allowed is 4".



Model 7202 Knee Space Drawers

Drawers have a fully enclosed shell with top cutouts for mounting [to countertops], zinc-plated flush pull handles, and full extension slides.

Single knee space drawer:

7202-18 = 18" Wide 7202-24 = 24" Wide

Double knee space drawer:

7202-36 = 36" Wide 7202-48 = 48" Wide

*Add drawer **ENDS** (attached to either or both sides of drawer housing) to accommodate Legs or fill additional space between cabinets. Available in any width up to 24"; specify width for *each* end (Ex.: 4" Left; 6" Right).

7202-30 = 24" Wide Single Drawer with (2) 3" Ends
 7202-30R = 24" Wide Single Drawer with 6" Right End

End Panels

As an option to H-Legs, End Panels are ideal for end-of-run applications, and especially where undercounter appliances are used. A built-in toe space gives the appearance of a finished cabinet side for a uniform appearance. Panels are 1" wide and available in two heights and two depths.

Standing Height 35-3/4" Sitting Height 28-1/2"

L-22 = 22-1/2" Deep LL-22 = 22-1/2" Deep
 L-29 = 29" Deep LL-29 = 29" Deep

Knee Space Panels

Expandable Knee Space Panels close pipe chase at rear of knee space areas. Panels come in four heights and two expansions: 22" to 36" wide; or 36" to 53" wide (Custom lengths or fixed panels are available upon request).


	<u>Standing Height</u>	<u>Sitting Height</u>
Full Height:	UN1000-01 (22"-36")	UN1000-05 (22"-36")
	UN1000-02 (36"-53")	UN1000-06 (36"-53")
Under 7200:	UN1000-03 (22"-36")	UN1000-07 (22"-36")
	UN1000-04 (36"-53")	UN1000-08 (36"-53")

Table Frames/Aprons/Legs

7200 Series – Table/Knee Space Components

Dimensions:

End View



Product Features:

- * Legs manufactured of 2" square steel tubing and include leveling feet.

Ordering Information:

- * Specify frame and apron widths in inches.
- * Optional Leg Boots (Model CMB-1) are available in black vinyl only.

Model 7210 Aprons

Aprons span the space between cabinets where full table frames are not necessary. Aprons are available in any width up to 66".

7210-** = Specify Width (Ex.: **7210-45** = 45" Wide)



Model 7210 Table Frames

Table Frames span the space between cabinets for stability and additional countertop support. Corner pockets are added at factory for Legs when required. Frames are available in any width up to 66".

7210-**(width) x 22 = 22-1/2" Deep


7210-**(width) x 28 = 28" Deep

(Ex.: **7210-35x22** = 35" Wide x 22-1/2" Deep)



Model 7228, 7236 Table Legs

H-legs and single legs can be used with table frames and knee space drawers (with ends). All Legs are available in two heights, and H-Legs are available in two depths.



Sitting Height 28-1/2"

Standing Height 35-3/4"

7228S = Single Leg

7236S = Single Leg

7228-22 = H-Leg 22-1/2"D


7236-22 = H-leg 22-1/2"D

7228-28 = H-Leg 28"D

7236-28 = H-Leg 28"D

*Use Leg Stretchers for added stability between H-Legs at independent tables.

Model LS-** (indicate length)



Fully assembled, your table will look like this. (Shown without countertop – see pg 41 for counter material options.)

Table Frames

Adjustable Height Tables - Sovella®


Product Features:

- * BASIC Lab Bench Frames are constructed of tubular steel. Tables are available with a wide variety of options and accessories. Any Table can be made mobile by adding casters.
- * Manual adjustment only.

Ordering Information:

- * Finish is standard ESD GRAY.
- * Contact factory for options and accessories.
- * Countertops are not included (see page 43).


BASIC TABLE



BASIC Adjustable Height Worktable Frames


Worktable frame size W x D inches	Worktable frame model	max load lbs
47.24 x 29.52	120 35 021	660
59.05 x 29.52	120 35 022	660
59.05 x 35.43	120 35 037	660
70.86 x 29.52	120 35 023	660
70.86 x 35.43	120 35 038	660
Height adjustment	25.59 to 35.43 inches	

SAMPLE CONFIGURATIONS:



Side-by-side Tables

Island Bench assembly with Basic Tables, countertop mounted upright shelving with brackets, power rails and 7300 series mobile cabinets



*See Page 47 for Upper Shelving

Table Frames

Adjustable Height Tables - Sovella®

Product Features:


- * CONCEPT Lab Bench Frames are ergonomically designed and available with a wide variety of options and accessories. Any Table can be made mobile by adding casters.
- * CONCEPT Tables are also available with a hand crank or motor adjustment.

- * Our 7300/7400 Line cabinets can be suspended from CONCEPT frames 29.53"D with manual adjustment only.


Ordering Information:

- * Finish is standard ESD GRAY.
- * Contact factory for options and accessories.
- * Countertops are not included (see page 43).

CONCEPT TABLE




CONCEPT Table with Optional casters




CONCEPT Adjustable Height Worktable Frames

Worktable frame size W x D inches	Standard table frame model	max load lbs	Handcrank height adjustment model	max load lbs	Motor height adjustment model	max load lbs
39.37 x 23.62	100 35 015	440	101 35 015	440	103 35 115	660
39.37 x 29.52	100 35 030	1100	101 35 030	440	103 35 130	660
47.24 x 23.62	100 35 016	440	101 35 016	440	103 35 116	660
47.24 x 29.52	100 35 031	1100	101 35 031	440	103 35 131	660
59.05 x 23.62	100 35 017	440	101 35 017	440	103 35 117	660
59.05 x 29.52	100 35 032	1100	101 35 032	440	103 35 132	660
70.86 x 23.62	100 35 018	440	101 35 018	440	103 35 118	660
70.86 x 29.52	100 35 033	1100	101 35 033	440	103 35 133	660
78.74 x 23.62	100 35 019	440	n/a		103 35 119	660
78.74 x 29.52	100 35 034	1100	n/a		103 35 134	660
Height adjustment	26.37--44.09 inches		26.57--40.35 inches		27.55-- 43.30 inches	



SAMPLE CONFIGURATIONS:

Mobile Concept Table with countertop, casters uprights, suspended cabinet, shelves, brackets, lower CPU holder and power rail




Concept Table with countertop, uprights, shelves, brackets, monitor arm, power rail

Wall Cabinets

7500, 7600, 7700, 7800 Series

Dimensions:

End View




Product Features:

- * 7500 Units have one shelf; 7600 and 7700 Units have two shelves; 7800 Units have three shelves.
- * Shelves are adjustable in 1/2" increments.
- * Glass is 1/4" thick.

Ordering Information:


- * Specify swing for doors on *single* units by substituting 'R' with 'R' or 'L' after series number (Units are shown with *left* swing doors). Ex. 7501**R**24 for right swing door.
- * For finished backs, add suffix 'FB' to model number.
- * Specify locks if required.

18" wide | 24" wide | 29" wide | 35" wide | 41" wide | 47" wide




Model 7*01 Hinged Solid Door Unit

24"H	7501-18	7501-24	-----	-----	-----	-----
30"H	7601-18	7601-24	-----	-----	-----	-----
36"H	7701-18	7701-24	-----	-----	-----	-----
48"H	7801-18	7801-24	-----	-----	-----	-----




Model 7*02 Hinged Solid 2-Door Unit

24"H	-----	-----	7502-29	7502-35	7502-41	7502-47
30"H	-----	-----	7602-29	7602-35	7602-41	7602-47
36"H	-----	-----	7702-29	7702-35	7702-41	7702-47
48"H	-----	-----	7802-29	7802-35	7802-41	7802-47



Model 7*03 Sliding Solid 2-Door Unit

24"H	-----	-----	7503-29	7503-35	7503-41	7503-47
30"H	-----	-----	7603-29	7603-35	7603-41	7603-47
36"H	-----	-----	7703-29	7703-35	7703-41	7703-47
48"H	-----	-----	7803-29	7803-35	7803-41	7803-47



Model 7*04 Hinged Framed Glass Door Unit


24"H	7504-18	7504-24	-----	-----	-----	-----
30"H	7604-18	7604-24	-----	-----	-----	-----
36"H	7704-18	7704-24	-----	-----	-----	-----
48"H	7804-18	7804-24	-----	-----	-----	-----

Wall Cabinets

7500, 7600, 7700, 7800 Series

Dimensions:

End View




Product Features:

- * 7500 Units have one shelf; 7600 and 7700 Units have two shelves; 7800 Units have three shelves.
- * Shelves are adjustable in 1/2" increments.
- * Glass is 1/4" thick.

Ordering Information:


- * For finished backs, add suffix 'FB' to model number.
- * Specify locks if required.



18" wide | 24" wide | 29" wide | 35" wide | 41" wide | 47" wide


Model 7*05 Hinged Framed Glass 2-Door Unit

	18" wide	24" wide	29" wide	35" wide	41" wide	47" wide
24"H	----	----	7505-29	7505-35	7505-41	7505-47
30"H	----	----	7605-29	7605-35	7605-41	7605-47
36"H	----	----	7705-29	7705-35	7705-41	7705-47
48"H	----	----	7805-29	7805-29	7805-41	7805-47



Model 7*06 Sliding Framed Glass 2-Door Unit

	18" wide	24" wide	29" wide	35" wide	41" wide	47" wide
24"H	----	----	7506-29	7506-35	7506-41	7506-47
30"H	----	----	7606-29	7606-35	7606-41	7606-47
36"H	----	----	7706-29	7706-35	7706-41	7706-47
48"H	----	----	7806-29	7806-35	7806-41	7806-47



Model 7*08 Open Unit

	18" wide	24" wide	29" wide	35" wide	41" wide	47" wide
24"H	7508-18	7508-24	7508-29	7508-35	7508-41	7508-47
30"H	7608-18	7608-24	7608-29	7608-35	7608-41	7608-47
36"H	7708-18	7708-24	7708-29	7708-35	7708-41	7708-47
48"H	7808-18	7808-24	7808-29	7808-35	7808-41	7808-47

Corner Wall Cabinets

7500, 7600, 7700, 7800 Series

Product Features:




- * 7500 (24"H) Units have one shelf; 7600 (30"H) and 7700 (36"H) units have two shelves; 7800 (48"H) units have three shelves.
- * Shelves are adjustable in 1/2" increments.
- * Glass is 1/4" thick.

Ordering Information:

- * Specify swing for doors by adding 'R' or 'L' after series number (Units shown with *left* swing doors). Ex. 7501RC for right swing corner unit.
- * Specify locks if required.


Dimensions:

End View



Model 7*01C
Hinged Solid
Door Unit

24"H	7501C
30"H	7601C
36"H	7701C
48"H	7801C



Model 7*04C
Hinged Framed
Glass Door Unit


24"H	7504C
30"H	7604C
36"H	7704C
48"H	7804C

Freestanding Tall Cabinets

7900, 7910 Series

Dimensions:

End View




Product Features:

- * Freestanders are available 18" or 24" deep.
- * All units include five adjustable shelves.
- * Shelves are adjustable in 1/2" increments.

Ordering Information:

- * Specify swing for doors on *single* units by substituting '-' with '**R**' or '**L**' after series number (Units are shown with *left* swing doors).
Ex. 7901**R**18 for right swing door unit, 18" wide).
- * Specify locks if required.



Model 7901

18"D Freestander includes five adjustable shelves behind hinged door.

7901-18 = 18" Wide


7901-24 = 24" Wide

Model 7911

24"D Freestander includes five adjustable shelves behind hinged door.

7911-18 = 18" Wide

7911-24 = 24" Wide



Model 7902

18"D Freestander includes five adjustable shelves behind hinged double doors.

7902-29 = 29" Wide

7902-41 = 41" Wide

7902-35 = 35" Wide

7902-47 = 47" Wide

Model 7912

24"D Freestander includes five adjustable shelves behind hinged double doors.

7912-29 = 29" Wide

7912-41 = 41" Wide

7912-35 = 35" Wide


7912-47 = 47" Wide

Freestanding Tall Cabinets

7900, 7910 Series

Dimensions:

End View




Product Features:

- * Freestanders are available 18" or 24" deep.
- * All units include five adjustable shelves.
- * Shelves are adjustable in 1/2" increments.
- * Glass is 1/4" thick.

Ordering Information:

- * Specify swing for doors on *single* units by substituting '-' with 'R' or 'L' after series number (Units are shown with *left* swing doors).
Ex. 7904R18 for right swing door unit, 18" wide.
- * Specify locks if required.



Model 7904

18"D Freestander includes five adjustable shelves behind hinged framed glass door.

7904-18 = 18" Wide


7904-24 = 24" Wide

Model 7914

24"D Freestander includes five adjustable shelves behind hinged framed glass door.

7914-18 = 18" Wide

7914-24 = 24" Wide



Model 7905

18"D Freestander includes five adjustable shelves behind hinged framed glass double doors.

7905-29 = 29" Wide

7905-41 = 41" Wide

7905-35 = 35" Wide

7905-47 = 47" Wide

Model 7915

24"D Freestander includes five adjustable shelves behind hinged framed glass double doors.

7915-29 = 29" Wide

7915-41 = 41" Wide

7915-35 = 35" Wide


7915-47 = 47" Wide

Freestanding Tall Cabinets

7900, 7910 Series

Dimensions:

End View




Product Features:

- * Freestanders are available 18" or 24" deep.
- * All units include five adjustable shelves.
- * Shelves are adjustable in 1/2" increments.
- * Glass is 1/4" thick.

Ordering Information:

- * Specify lock if required.



Model 7906

18"D Freestander include five adjustable shelves behind sliding framed glass doors.

7906-29 = 29" Wide

7906-41 = 41" Wide

7906-35 = 35" Wide

7906-47 = 47" Wide

Model 7916


24"D Freestander includes five adjustable shelves behind sliding framed glass doors.

7916-29 = 29" Wide

7916-41 = 41" Wide

7916-35 = 35" Wide

7916-47 = 47" Wide



Model 7908

18"D Open freestander includes five adjustable shelves.

7908-18 = 18" Wide

7908-35 = 35" Wide

7908-24 = 24" Wide

7908-41 = 41" Wide

7908-29 = 29" Wide

7908-47 = 47" Wide

Model 7918

24"D Open freestander includes five adjustable shelves..

7918-18 = 18" Wide

7918-35 = 35" Wide

7918-24 = 24" Wide

7918-41 = 41" Wide

7918-29 = 29" Wide

7918-47 = 47" Wide

Specialty Cabinets


Acid Storage, Narcotic Storage

Product Features:

- * Cabinets feature the same construction as standard cabinets with additional features as specified.
- * Shelves are adjustable in 1/2" increments.

Ordering Information:

- * For finished backs, add suffix 'FB' to model number.
- * ADA height unit also available, 32-1/8"H



Model 7000

Acid/Corrosive Resistant Vented Safety Cabinet **Model 7000**

The 7000 cabinet safely stores acids and solvents within a welded steel unit lined with Wesliner material for maximum corrosion resistance in the event of spills or accidents. Each of the two hinged doors is vented top and bottom to prevent any build up of noxious fumes.

The Wesliner lining used is a monolithic, natural white non-asbestos fiber/cement panel with a smooth surface on both sides. It is a strong and rugged material designed to offer maximum resistance to heat, chemicals, impact, and moisture.

Unit is standard 35-3/4"H x 22-1/2"D and includes one full width white Phenolic resin adjustable shelf. Standard "ACID" label is also included (red letters on silver).

7000-29 = 29" Wide

7000-35 = 35" Wide

7000-41 = 41" Wide


7000-47 = 47" Wide

Narcotic Locker Cabinet

Model 7601-24NC

Model 7601-24NC narcotic locker features vermin-free double wall construction, two adjustable shelves, dual locking double door assembly, and lighted alarm for maximum security of floor stock class narcotic medications. Unit can be wall or counter mounted.

The dual locking double door assembly consists of an 18 gauge double pan welded door and a full height 1/8" plate inside door suspended on a piano hinge welded to the cabinet to eliminate removal of the hinge pins. Both doors are designed with interchangeable removable core locks* and are complimented with electronic circuitry. A warning light will automatically be activated if the inside door is left unlocked (the light can be remotely wired to a Director's office to alert entry to the locker). An optional audible alarm is available upon request.



Model 7601-24NC

Model **7601-24NC** = 24"W x 30"H x 13"D (standard)

*Core locks have 50 core locking combinations. Master keying is available.

Wall Cabinet/Freestander Accessories

Sloping Tops


Product Features:

- * Sloping Tops include end caps on both sides.
- * Soffit panels available upon request.



Ordering Information:

- * Specify special length sloping tops in inches (Ex. 7613-44 for 44" top length).

End Views:




Model 7613



Model 7618

Depth	18" wide	24" wide	29" wide	35" wide	41" wide	47" wide
↓						
13"	7613-18	7613-24	7613-29	7613-35	7613-41	7613-47
18"	7618-18	7618-24	7618-29	7618-35	7618-41	7618-47
24"	7624-18	7624-24	7624-29	7624-35	7624-41	7624-47




Model 7624

Corner Sloping Tops

Corner tops used where wall cabinets meet at a 90° corner.


Depth	Model
↓	
13"	7613C
18"	7618C
24"	7624C



Corner Cabinet Sloping Top

Corner top for use with corner wall cabinet series 7500C, 7600C, 7700C, or 7800C only.

Model 7613CC = 13" Deep x 18" Front



Cabinet Accessories


Fillers

Product Features:

- * Finish is the same as cabinets.
- * Concealed clips hold fillers in place.

Ordering Information:

- * Specify filler height and width in inches
(Ex. FBF36-6½ = Front Base Filler 35-3/4”H x 6-1/2” Wide).



Model FBF

Front Base Filler with toe space is used between adjacent base cabinets, or cabinets adjacent to a wall or end sink cabinet.

FBF36 = Standing Height 35-3/4”H


FBF28 = Sitting Height 28-1/2”H

Model SBF

Side [End] Base Filler is used to close the chase between the base cabinet and the wall, or back-to-back cabinets at islands/peninsulas.

SBF36 = Standing Height 35-3/4”H


SBF28 = Sitting Height 28-1/2”H



Model CTF-4

Corner Toe Filler is used to close the toe space opening where base cabinets meet at a 90° corner.

CTF-4 = 4” High



Model WCF

Wall Cabinet Filler is used between adjacent wall cabinets, or cabinets adjacent to a wall. Fillers are standard 13” Deep.

WCF24 = 24” High (for 7500 Series cabinets)

WCF30 = 30” High (for 7600 Series cabinets)

WCF36 = 36” High (for 7700 Series cabinets)

WCF48 = 48” High (for 7800 Series cabinets)

Cabinet Accessories


Risers, Cabinet Shelves, and Lazy Susans

Product Features:

- * Riser Panels are sold in pairs.
- * Cabinet shelves over 36" have reinforcing channels spot-welded to underside.

Ordering Information:

- * All 199 shelves are finished in *Champagne*, unless otherwise specified when ordering.
- * Lazy Susan shelf colors are as noted.



Model 260

Riser panels are sold as a pair, and attach to either side of a wall unit. 44" High x 13-7/8" Deep (overall). Designed for use with 13" Deep wall units only.




Model 199 - Cabinet Shelves

Adjustable shelves include four mounting clips each. Specify model number of unit to receive shelf after the shelf model number (Ex. **199-54 x 7002-18** = Base unit shelf for 18" Wide 7002 model).

Model 199-54 = Base Unit Shelf

Model 199-55 = Wall Unit Shelf

Model 199-79 = Freestander Shelf



Lazy Susan

Each Lazy Susan set includes two 24" Plastic Polymer "Dependently-Rotating Shelves" (both shelves move at the same time) and telescoping chrome-plated steel center post.

*For use with Corner Base Unit Models 7171/7471 only.

Model REV-6012-24-**A** = Lazy Susan in **Almond**

Model REV-6012-24-**W** = Lazy Susan in **White**

Laboratory Accessories

Miscellaneous Supports

Product Features:

- * Pipe frame assemblies have pre-galvanized finish unless otherwise noted.
- * Angles are finished in Champagne.


Ordering Information:

- * Specify height and width of Pipe support frames in inches (Ex. **UNI36-8** = Pipe frame 35-3/4"H x 8"W)
 - * Specify custom lengths of wall angles in inches.
-




Model UNI

Pipe Frame Support Assemblies are used to support plumbing lines and worktop overhangs. Assembly includes (2) vertical channel struts with 14 gage bent "U" bracket at top and bottom. Assemblies are sized to fit the service chase behind base cabinets at wall benches, islands and peninsulas. Pipe hangers are not included (supplied by plumber).



Model WA

Wall angles provide support for countertops where legs or end panels are not needed or cannot be used. Wall angles are 2" x 2" x 18" long, 14 gauge. Special lengths are available upon request.



Model CWA

Corner wall angles are used to help support countertops behind corner base units or 90° dead space areas. Corner wall angles are 2" x 2" x 18" long, 14 gauge. Special lengths are available upon request.

Laboratory Accessories


Wall Shelving

Product Features:

- * Plastic Laminate and Resin shelves are finished on all edges, top, and bottom.
- * Stainless steel shelves have bottom edge returns for safety. 18ga Type 304 with #4 finish.

Ordering Information:

- * Brackets and standards are ordered by each piece (not sold in pairs).
- * Consult factory for availability of colors, sizes, or other bracket/standard types not listed.



Model SH*

Wall shelves are available in any size required. Specify size in inches.

SH1 = Plastic Laminate


SH2 = Epoxy Resin

SH3 = Stainless Steel

SH4 = Lab Grade Plastic Laminate

SH8 = Phenolic Resin


Example: Model **SH2-3612** = Epoxy shelf 36"W x 12"D



Model 85ANO

Heavy Duty Double slotted standard for use with Model 185 bracket. Standard is 1-1/4" wide x 1/2" high with 1" slot adjustment. Available lengths: 24", 36", 48", 60", 72", 84", and 96". 16 gage, Anochrome finish standard.


Example: Model **85ANO-36** = 36" Long Standard



Model 185ANO

Heavy duty double bracket for use with Model 85ANO standard. Anochrome finish standard. Available depths: 8", 10", 12", 14", 16", 18", 20", 22", and 24". 16 gage for brackets up to 12", 14 gage for brackets 14" to 24", Anochrome finish standard.

Example: Model **185ANO-12** = 12"D Bracket



Model 106ANO

Shelf Rest is used for joining two shelves at one bracket. Provides positive mounting of shelves to brackets. Use with 185ANO brackets only.



Model PB

Book-end style shelf brackets are cold rolled steel and painted to match cabinetry. Can be used with 85ANO standards. Ends are 6"H, depth made to match shelf. One each left and right bracket are required for each shelf.

Example: Model **PB-12** = 12"D Single Bracket

Laboratory Accessories

Island/Peninsula Upright Shelving



Product Features:

- * Double sided/double slotted tubular steel upright posts can be mounted to the floor or the countertop.
- * Various steel brackets are available in a variety of styles.
- * Any shelf material can be used (painted or stainless steel, epoxy, laminate, etc.)



Ordering Information:

- * Consult factory for other configurations, post heights, bracket styles, etc.
- * Steel components can be painted to match the casework.

SAMPLE CONFIGURATIONS



'H'-Style Double sided island shelving assembly with book end brackets, center shelves, countertop mounted posts with foot plates.



Laboratory Accessories

Reagent Rack Assemblies, Service Chases


Product Features:

- * Steel uprights and service chases have a removable access panel for easy plumbing access.
- * Custom assemblies can be designed upon request to meet laboratory requirements.

Ordering Information:

- * Specify dimensions required for uprights, fascias, and service chases in inches.
- * Specify fixed or adjustable intermediate reagent shelves.
- * Specify shelf size and material.

TYPICAL CENTER BENCH or PENINSULA DOUBLE TIER REAGENT ASSEMBLY:



Rack Components

- RRU** Single tier END upright (4"D standard)
- RRUC** Single tier CENTER upright (4"D standard)
- RRDU** Double tier END upright (4"D standard)
- RRDUC** Double tier CENTER upright (4"D standard)
- RRFP** Reagent rack fascia panel mounts underneath shelf between uprights, for added strength and rigidity (3"H standard)
- RRFC** Reagent rack fixture is mounted underneath shelf for pendant-hung [under-mounted] service fixtures
- VSCE** Vertical Service Chase Enclosure (not shown) mounts above countertop and extends to or above ceiling; includes removable access panel on one side

Shelves

- CT-1RS** Plastic Laminate Reagent shelf 1-1/4" thick for top shelf; 1-1/4" or 3/4" thick for intermediate shelves
- CT-4RS** Lab Grade Plastic Laminate Shelf 1-1/4" thick for top shelf; 1-1/4" or 3/4" thick for Intermediate shelves
- CT-2RS** Epoxy resin shelf 1" thick
- CT-8RS** Phenolic resin shelf 1" thick

Laboratory Equipment

Fume Hood Base Cabinets


Base cabinets are available to fit under standard laboratory work surfaces and counters, with safety features suitable for hazards as specified.

Standard cabinets follow the same design and construction as 7000 series cabinets.

Acid cabinets follow the same design and construction as 7000 series cabinets but include an interior lining of Wesliner material and front louvers for venting. "ACID" label is included. (See page 39 for more info)


Flammable Storage cabinets feature double wall construction and rear vents. Available with Manual or Self-close (Sliding) doors. OSHA and NFPA compliant. (Units are finished in SAFETY YELLOW unless otherwise specified)

Standard




Model	Dim. W x H
7119-29FH	29" x 35-3/4"
7119-35FH	35" x 35-3/4"
7119-47FH	47" x 35-3/4"
7119-58FH	58" x 35-3/4"

Acid Storage



Model	Dim. W x H
7000-29FH	29" x 35-3/4"
7000-35FH	35" x 35-3/4"
7000-47FH	47" x 35-3/4"
7000-58FH	58" x 35-3/4"

Flammable Storage



Model	Dim. W x H x D
L124 Manual	35" x 35-9/16" x 22"
L224 Self-Close	35" x 35-5/8" x 22"
L136 Manual	47" x 35-9/16" x 22"
L236 Self-Close	47" x 35-5/8" x 22"
L144 Manual	59" x 35-9/16" x 22"
L244 Self-Close	59" x 35-5/8" x 22"

Vent Kits

Model VK-A

Vent Kit for 7000 series Acid cabinets. PVC Vent extends from the cabinet to the worktop so fumes can be exhausted by the hood. Kit includes elbow and two threaded connections.

Model VK-F

Vent Kit for Flammable storage cabinet. Flexible plastic hose extends from the cabinet directly into ductwork, bypassing the fume hood entirely so dangerous fumes are exhausted directly into the ductwork. Kit includes 10' hose and adapters.

*Other materials available upon request.

